

## Table of contents

1. Reagents and equipment .....	2
2. Passivation protocol .....	2
2.1. Choosing the right passivation .....	2
2.2. Passivation with PLL-PEG .....	3
2.3. Passivation with PEG-SVA .....	3
2.3.1. PLL or laminin-based passivation.....	4
2.4. APTES based passivation .....	4
3. Use of the PLPP photoinitiator .....	4
4. Photopatterning with PRIMO .....	4
5. Cell seeding .....	5
6. Related publications .....	5



This protocol is a guideline and you might need to optimize it for each different application.

## 1. Reagents and equipment

### Equipment

- Plasma cleaner (optional depending on the passivation)
- PRIMO® setup
- PDMS stencil

### Reagents

- PLPP 1X
- PBS 1X
- Protein solution of interest (typical concentration: 10-100  $\mu\text{g}\cdot\text{ml}^{-1}$ )

Some specific reagents are needed depending on the passivation. Please refer to the paragraph [2.1](#) for advice to choose your passivation.

#### PLL-g-PEG passivation:

- PLL-g-PEG (*Alvéole*, PLL(20)-g[3.5]-PEG(2) 0,1mg.ml<sup>-1</sup> in PBS 1x, storage : -20°C, 10x solution)

#### PEG-SVA passivation:

- mPEG-SVA (mPEG-Succinimidyl Valerate, MW 5 000, *Laysan bio*)
- HEPES 8,3 < pH < 8,5 at 0,1 M
- Parafilm®
- H<sub>2</sub>O mQ
- EtOH (70% or more)
- Surfactant Mix (only for soft PDMS, *Alvéole*)

#### PLL + PEG-SVA passivation:

- Polylysine (PLL, *Merck/Sigma-Aldrich*, ref: P8920)

#### Laminin + PEG-SVA passivation :

- Laminin

#### APTES + PEG-SVA passivation :

- 3 aminopropyltriethoxysilane (APTES, *Sigma-Aldrich* ref A3648).

## 2. Passivation protocol

### 2.1. Choosing the right passivation

#### PLL-g-PEG

The easiest way to passivate. **Works only on glass substrates and EM grids.** Using a plasma cleaner (or a super clean coverslip) is mandatory. The passivation is less strong than the following ones. Suitable for 24 hours to 48 hours experiments.

### PEG-SVA based passivation

These kinds of passivation are stronger and are more stable in time. Your cells will be confined on the patterned area for days.

#### PLL + PEG-SVA

The PLL + PEG-SVA passivation by default. Easier and cheaper than the other PEG-SVA passivation, it works for most applications all substrates.

#### Laminin + PEG-SVA

This passivation is more suited for patterning cardiomyocytes, induced pluripotent derived stem cells (IPS) and embryonic stem cells (ES).

#### APTES + PEG-SVA

Recommended for microtubules and purified protein patterning. It also works for classic cell lines and could also be an alternative for patterning cardiomyocytes, IPS and ES.

## 2.2. Passivation with PLL-PEG

*This passivation works only on glass substrates and EM grids*

Prepare the PDMS stencil by removing areas corresponding to the well.

- Introduce your substrate in the plasma cleaner. Remove any lid that could prevent the plasma cleaning process.
- Plasma air 1 minute.
- Place the PDMS stencil on the glass surface.
- Add 10-15  $\mu\text{L}$  of PLL-g-PEG at 100  $\mu\text{g}\cdot\text{ml}^{-1}$  in the wells.


*Important: from this point, never let the surface dry or dewet.*

- Incubate for 1 hour at room temperature. Cover your substrate or put it in a humid chamber to avoid evaporation. Be careful of drying.
- Rinse the wells 4 to 5 times with PBS 1x.
- Leave 10-15  $\mu\text{L}$  of PBS 1x in each well.

## 2.3. Passivation with PEG-SVA

The strategy for passivation with PEG-SVA is the following:

- 1) Coating of the substrate with PLL, APTES or laminin.
- 2) Coupling of the amines of the PLL, APTES or laminin with the SVA ester.



mPEG-SVA is a PEG (5kDa) coupled with a Succinimidyl Valerate (SVA) ester. This ester is used for PEGylation of primary amine.

### 2.3.1. PLL or laminin-based passivation

Prepare the PDMS stencil by removing areas corresponding to the well.

- Plasma clean your substrate for 1 min. Optional, but highly recommended.
- Put the PDMS stencil on the surface of your substrate.
- Add 15  $\mu\text{l}$  of PLL at 100  $\mu\text{g}\cdot\text{ml}^{-1}$  diluted in mQ water for 30 minutes or laminin at 100  $\mu\text{g}\cdot\text{ml}^{-1}$  in PBS 1 X in the wells.

*Important: from this point, never let the surface dry or dewet.*

- Incubate 1 hour at room temperature. Cover your substrate or put it in a humid chamber to avoid evaporation.
- Rinse the wells 3 to 4 times with HEPES.
- Make a solution of PEG-SVA in 0,1 M HEPES with  $8,3 < \text{pH} < 8,5$  at 70  $\text{mg}\cdot\text{ml}^{-1}$ . If the passivation is not strong enough, you can increase the concentration up to 100  $\text{mg}\cdot\text{ml}^{-1}$ .

**As the half-life of the SVA ester is 10 minutes at pH 8,5, the solution must be prepared just before use.**

- Add 10-15  $\mu\text{l}$  of the PEG-SVA solution. Incubate 1 hour at room temperature, in humid chamber to avoid evaporation. Be carefull of drying.
- Rinse the wells 4 to 5 times with PBS 1x.
- Leave 10-15  $\mu\text{L}$  of PBS 1 X in each well.

### 2.3.2. APTES-based passivation

- Plasma clean your substrate for 1 min.
- Transfer your substrate in a plastic box (plastic petri dish for example).

*Silanization must be performed under a chemical hood*

- Detach the cap of a 1,5 mL Eppendorf 1,5 mL tube and put it in the plastic box near to the substrate.
- Fill the cap with 200  $\mu\text{l}$  of APTES.
- Close the box and seal it with Parafilm<sup>®</sup>. Wait 1 hour.
- Put the PDMS stencil on the surface of your substrate.
- Make a solution of PEG-SVA in 0,1M HEPES with  $8,3 < \text{pH} < 8,5$  at 70  $\text{mg}\cdot\text{ml}^{-1}$ . If the passivation is not strong enough, you can increase the concentration up to 100  $\text{mg}\cdot\text{ml}^{-1}$ .

**As the half-life of the SVA ester is 10 minutes at pH 8,5, the solution must be prepared just before use.**

- Incubate with 15  $\mu\text{l}$  of the PEG-SVA solution for one hour. Cover your substrate or put it in a humid chamber to avoid evaporation.

*Important: from this point, never let the surface dry or dewet.*

- Rinse the wells 4 to 5 times with PBS 1x.
- Leave 10-15  $\mu\text{L}$  of PBS 1x in each well.

## 3. Use of the PLPP photoinitiator

- Remove the solution from the well without drying the surface.
- Fill the wells with 5-15 ml of PLPP.

## 4. Photopatterning with PRIMO

- Place your substrate on the microscope.
- Check that the PRIMO UV filter is on the optical pathway.
- Turn on the PRIMO set up.



- Load and lock your pattern on Leonardo software. The typical UV doses are listed below. It might depend on the density of PLPP Gel, substrate, etc... **dose adjustment might be needed.**

	PRIMO (laser 375 nm)	PRIMO 2 (LED 365 nm)
PLL-g-PEG	600 mJ.mm <sup>-2</sup>	300 mJ.mm <sup>-2</sup>
PEG-SVA	1200 mJ.mm <sup>-2</sup>	600 mJ.mm <sup>-2</sup>

- Focus on the surface you want your pattern to be. This step is crucial to have nice and sharp pattern. You can use your substrate imperfection like dusts to focus.
- Launch patterning sequence.
- Rinse the wells 3 times with 20 µL of PBS 1x.
- Remove the solution from the wells but leave a few microliters so the surface will not dry out.
- Incubate 5 minutes with a solution of 20 µg/ml of protein in its dedicated buffer. You can increase the concentration to 100 µg.ml<sup>-1</sup> and the incubation time to 15 min.

**Longer incubation time will lead to higher protein concentration on patterns but also higher background.**

- Rinse the well 4-5 times with 20 µL of PBS 1 X
- Leave 10-15 µL of PBS 1 X in the well.

## 5. Cell seeding

This step is really cell type dependent and will need some adjustment.

- Rinse your sample with media without dewetting the substrate.
- Add 1.10<sup>4</sup> to 1.10<sup>5</sup> cells per cm<sup>2</sup>
- Carefully put your sample in the incubator and let them adhere for 15 min to 1 hour (or more) depending on your cell type
- Rinse gently with cell culture media to remove the non-adherent cells.
- If cells are attaching outside the patterns, wait one more hour and pipet over the cells until the one outside of the pattern detach.

## 6. Related publications

### LIMAP technology:

- *Multiprotein Printing by Light-Induced Molecular Adsorption*, Strale et al., *Advanced materials*, 2015.

### Classical proteins photopatterning (glass with PLL-g-PEG):

- *Altered microtubule dynamics and vesicular transport in mouse and human MeCP2- deficient astrocytes*, Delépine et al., *Hum Mol Genet.*, 2016.
- *Collagen assembly and turnover imaged with a CRISPR-Cas9 engineered Dendra2 tag*, Pickard et al., *bioRxiv*, 2018.
- *The cytoskeleton as a smart composite material: A unified pathway linking microtubules, myosin-II filaments and integrin adhesions*, Rafiq et al., *bioRxiv*, 2018.

### Patterning on soft PDMS:

- *Extracellular matrix stiffness regulates force transmission pathways in multicellular ensembles of human airway smooth muscle cells*, Polio et al., *bioRxiv*, 2018.

Proteins photopatterning on 3D structures:

- *A new approach to design artificial 3D micro-niches with combined chemical, topographical and rheological cues, Stoecklin et al., Adv biosys, 2018.*
- *Controlling the topography and biochemistry of cell culture substrates with PRIMO® photopatterning system, Strale et al., Nature Methods, 2018 (Application note).*

