

# PRIMO PROTOCOL: HYDROGEL PHOTOPOLYMERIZATION

Rev 2 – last update 15.01.2018

This protocol is a guideline and you might need to optimize it for each different application.

## Protocol for the use of 4 ARM PEG-ACRL, MW 10,000 hydrogel (Laysan Bio, Inc.)

### Equipment

- Inverted microscope
- Motorized XY-stage
- Camera
- PRIMO®
- Control software LEONARDO
- Computer
- Micropipettes and adequate tips

### Reagents

- 4 ARL PEG-ACRL, MW 10,000 hydrogel
- PLPP 1x solution
- Deionized water
- PBS 1x solution pH = 7,4

## 1. Hydrogel solution preparation

4 ARL PEG-ACRL, MW 10,000 hydrogel

- Dissolve hydrogel at 10% in H<sub>2</sub>O.
- Add 20 to 50 % of PLPP to the hydrogel solution.
- Wash 5 times with PBS.

## 2. Illumination

The height of the hydrogel will depend on the height of illuminated liquid. To control it is possible to place the hydrogel solution in a chamber of controlled height.

- Place the solution in the PDMS chamber.
- Focus in your well.
- Illuminate the hydrogel with your UV pattern (typical UV dose: 300-900mJ.mm<sup>-2</sup>).

***If working with NIS plugin insolate for 10 to 30s***

- Rinse 5 times with PBS.

You can control the reticulation of the gel by changing the hydrogel concentration, the PLPP percentage and the UV dose.

# PRIMO PROTOCOL: HYDROGEL PHOTOPOLYMERIZATION

Rev 2 – last update 15.01.2018

Protocol for the use of Collagen I – Methacrylate (Photocol kit, Methacrylated Type I Collagen for crosslinkable hydrogels, cat n°5201-1EA, lot n°7519, Advanced Biomatrix).

## Equipment

- Inverted microscope
- Motorized XY-stage
- Camera
- PRIMO®
- Control software LEONARDO
- Computer
- Micropipettes and adequate tips

## Reagents

- Collagen I – Methacrylate (Advanced Biomatrix)
- PLPP 1x solution
- Acetic acid
- PBS 1x solution pH = 7,4

## 3. Hydrogel solution preparation

Collagen I – Methacrylate (Advanced Biomatrix)

- Dissolve hydrogel at 8mg/mL in acetic acid.
- Add 50 % of PLPP to the hydrogel solution.
- Micropipettes and adequate tips
- Wash 5 times with PBS.

## 4. Illumination

The height of the hydrogel will depend on the height of illuminated liquid. To control it is possible to place the hydrogel solution in a chamber of controlled height.

- Place the solution in the PDMS chamber.
- Focus in your well.
- Illuminate the hydrogel with your UV pattern typical UV dose: 300-600mJ.mm<sup>-2</sup>).
- Rinse 5 times with PBS.

You can control the reticulation of the gel by changing the hydrogel concentration, the PLPP percentage and the illumination time. **!/ Always work on ice.**

# PRIMO PROTOCOL: HYDROGEL PHOTOPOLYMERIZATION

Rev 2 – last update 15.01.2018

Protocol for the use of 4 ARM PEG-ACRL, MW 10,000 hydrogel (Laysan Bio, Inc.) AND Collagen I – Methacrylate (Photocol kit, Methacrylated Type I Collagen for crosslinkable hydrogels, cat n°5201-1EA, lot n°7519, Advanced Biomatrix).

## Equipment

- Inverted microscope
- Motorized XY-stage
- Camera
- PRIMO®
- Control software LEONARDO
- Computer
- Micropipettes and adequate tips

## Reagents

- 4 ARL PEG-ACRL, MW 10,000 hydrogel
- Collagen I – Methacrylate
- PLPP 1x solution
- Deionized water
- Acetic acid
- PBS 1x solution pH = 7,4

## 5. Hydrogel solution preparation

4 ARL PEG-ACRL, MW 10,000 hydrogel

- Dissolve hydrogel at 10% in H<sub>2</sub>O.
- Add 20 to 50 % of PLPP to the hydrogel solution.
- Micropipettes and adequate tips
- Wash 5 times with PBS.

Collagen I – Methacrylate (Advanced Biomatrix)

- Dissolve hydrogel at 8mg/mL in acetic acid.
- Add 50 % of PLPP to the hydrogel solution.
- Micropipettes and adequate tips
- Wash 5 times with PBS.

For the mix of the both hydrogels:

- Optimal concentration: 80% Collagen I – Methacrylate (8mg/mL) + 20% 4arm-PEG-Acrylate (10%)
- Add 50% of PLPP to the hydrogel mix

## 6. Illumination

The height of the hydrogel will depend on the height of illuminated liquid. To control it is possible to place the hydrogel solution in a chamber of controlled height.

- Place the solution in the PDMS chamber.

## PRIMO PROTOCOL: HYDROGEL PHOTOPOLYMERIZATION

---

Rev 2 – last update 15.01.2018

- Focus in your well.
- Illuminate the hydrogel with your UV pattern (typical UV dose: 300-600mJ.mm<sup>-2</sup>).
- Rinse 5 times with PBS.

You can control the reticulation of the gel by changing the hydrogel concentration, the PLPP percentage and the illumination time. **!\ Always work on ice and prepare the mix extemporaneously.**

